



UmanoMedical

COMFORT SERIES

comfortnest greenXi™




Pressure Redistribution • Comfort • Durability

KEY BENEFITS

- 3 layer / 3 zone construction with varying density and IFD ratings for optimal pressure redistribution in the higher risk areas
- One-piece cover integrates all mattress components
- Viscoelastic memory gel foam topper for optimal immersion and comfort
- Entire foam core is constructed of 100% medical grade non-allergenic foam
- All components are anti-microbial and assist in prohibiting growth of bacteria and fungi
- Clam shell design cover features double-folded seams and 4 reinforced handles
- 1000 lb weight capacity

PRODUCT FEATURES



TECHNICAL SPECIFICATIONS

Dimensions	Center mattress: 35" W x 80" or 84" L x 6" H (88.9 cm W x 203.2 cm or 213.4 cm L x 15.2 cm H) Two lateral bolsters: 6" W x 80" or 84" L x 6" H (15.2 cm W x 203.2 cm or 213.4 cm L x 15.2 cm H) Foot end extender: 35" W x 4" L x 6" H (88.9 cm W x 10.2 cm L x 15.2 cm H)
Weight Capacity	1000 lb (454 Kg)
Flammability	Meets Federal Standard 16 CFR 1632, Cal. TB 117 Optional barrier to meet Cal. TB 129 and Boston Fire IX-11
Guarantee	Non prorated 5-year guarantee

Umano Medical reserves the right to change specifications without notice.



Umano Medical Inc.
230, boul. Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll Free: 1 844 409-4030
information@umanomedical.com
umanomedical.com