



umano medical

COMFORT SERIES

comfortnest yellow™




Comfort • Durability • Performance

KEY BENEFITS

- 1 layer / 1 zone reversible design (head to toe)
- Entire foam core is constructed of 100% medical grade non-allergenic foam
- All components are anti-microbial and assist in prohibiting growth of bacteria and fungi
- Clam shell design cover features double-folded seams and 4 reinforced handles
- 450 lb weight capacity
- comfortnest yellow+ model offers High Resilience (HR) foam and a weight capacity of 950 lb

PRODUCT FEATURES



TECHNICAL SPECIFICATIONS

Dimensions	35" W x 80" or 84" L x 6" H (88.9 cm W x 203.2 cm or 213.4 cm L x 15.2 cm H) Custom sizes available
Specialty Model	Behavioral & Mental Health model available: comfortnest yellow ^{MH} Fully enclosed cover, no handles, no flaps, no zippers.
Weight Capacity	450 lb (204 Kg) 950 lb (432 kg) on the comfortnest yellow+ model
Flammability	Meets Federal Standard 16 CFR 1632, Cal. TB 117 Optional barrier to meet Cal. TB 129 and Boston Fire IX-11
Guarantee	Non prorated 5-year guarantee on the comfortnest yellow and comfortnest yellow+ Non prorated 3-year guarantee on the comfortnest yellow ^{MH}

Umano Medical reserves the right to change specifications without notice.



Umano Medical Inc.
230, boul. Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll Free: 1 844 409-4030
information@umanomedical.com
umanomedical.com