



umano medical

THERAPEUTIC SERIES

covernest gray™




## Microclimate Management

### KEY BENEFITS


- Vapor permeable triple layered top sheet can be placed over any non-powered pressure redistribution surface (foam or gel surface) to add low air loss therapy.
- The low air loss feature addresses maceration and heat build up to help manage the patient's skin microclimate environment.

### Therapies

- On-Demand Low Air Loss



PRODUCT FEATURES



TECHNICAL SPECIFICATIONS

Available Sizes	Width: 36" (91.4 cm)
	Length: 80" and 84" (203 cm and 213 cm)
	Height: Fits 7" to 10" surfaces (17.8 cm to 25.4 cm)
Weight Capacity	as per weight capacity of surface used with LAL Cover
Control Unit	Dimensions 13.8" x 11" x 6.5" (35 cm x 28 cm x 16.4 cm)
	Weight 10.2 lb (4.6 kg)
	Power Cord 14' hospital-grade detachable
	Input Voltage 90 264 V, 50/60 Hz, 1A/0.5A, Universal
	ETL Listed Conforms to UL STD 60601-1 with respect to Electrical Shock, Fire and Mechanical Hazards
Certified to	CAN/CSA STD C22.2 No. 601.1
Flammability	Meets flammability standard CAL 117

Specifications can change without notice.



**Umano Medical Inc.**  
230, boul. Nilus-Leclerc  
L'Islet QC G0R 2C0 CANADA  
Toll Free: 1 (844) 409-4030  
information@umanomedical.com  
umanomedical.com

Umano Medical is distributor of the Therapeutic Series manufactured by KAP MEDICAL