



umano medical

BARIATRIC & MED-SURG BED
INFINITE™ WIDTH ADJUSTMENT
1060 LB CAPACITY



Unique design approach maximizes patient safety and dignity

Infinite™ Technology and built-in scale allow to reduce patient handling while improving staff efficiency and patient comfort

1060 lb capacity and 11" low height allow for an increased versatility throughout environments

VERSATILITY • LOW HEIGHT • INFINITE™ WIDTH ADJUSTMENT



MEDICAL-SURGICAL
Low height and fully equipped design maximize performance in med-surg units



BARIATRIC
Infinite™ width expansion and higher weight capacity suit the needs of bariatric care settings



PALLIATIVE CARE
Certification for use by two people, one patient and one bystander, ensures comfort in palliative care environments



MAXIMIZED SAFETY
Infinite™ width maintains a gap of under 60 mm between the head board and siderails
In compliance with standard IEC 60601-2-52



umano medical

BARIATRIC & MED-SURG BED

INFINITE™ WIDTH ADJUSTMENT

1060 LB CAPACITY

BED FEATURES

- 4 DC actuators
- Battery back-up
- Infinite™ width adjustment
- Compatibility with mattresses within 48" width
- 11" low height with 5" dual casters
- Built-in-bed scale system
- Programmable lowest height
- Programmable Auto-Contour mode
- Built-in bed exit detection system - 3 levels of sensitivity
- Remote bed exit alarm connectivity
- Smart screen, nurse controls and patient control lockout on footboard
- Nurse and patient controls integrated into siderails
- Trendelenburg angle displayed on footboard
- One hand release siderails
- Integrated bed extender (35"x 4")
- CPR release handle
- Lifiable and flat deck covers
- One step Trendelenburg
- One step vascular foot position
- Integrated pump holder
- Boostless™ Backrest System
- Backrest angle indicator
- Nightlights
- Bed exit side view lights
- Patient restraint loops
- 120-volt / 5A insulated outlets (at foot end)
- Communication port (for upgradable bed software)
- One-touch button groups
 - 30-degree head of bed
 - One-button chair position
 - Flat bed position

OPTIONAL FEATURES

- 6" dual casters
- Downward motion interrupt
- Bed status monitoring
- Smart Nightlights™
- Optional patient surface length (80" or 84")
- 120-volt/5A insulated outlets with USB mobile device power outlets
- Higher siderails (for 8" mattress)
- Foot end back-up control on the frame
- Different styling options available
- Roller type wall protection bumpers
- Additional nurse call port
- One-touch button groups
 - Backlit patient controls
 - Nurse call

Accessories

- Patient helper system (300 lb capacity)
- Upright oxygen bottle holder
- IV poles (permanent & removable)
- 4" bed extender cushion
- Removable patient controls
- Transport cart

Specifications

ook«SNOW» ALL

Product model	200-0000
Overall length (80" deck length)	90" / 228.6 cm
Overall length (84" deck length)	94" / 239 cm
Extended length	Adds 4" / 10.1 cm
Overall width	
Siderails down - 35" surface	39.75" / 101 cm
Siderails up - 35" surface	41.5" / 105.4 cm
Siderails down - 47" surface	51.5" / 130.8 cm
Siderails up - 47" surface	52.75" / 134 cm
Safe Working Load	1060 lb / 481 kg
Height range (to litter top)	
Low (with 5" / 12.7 cm dual casters)	11" / 27.94 cm
High (with 5" / 12.7 cm dual casters)	30" / 76.2 cm
Low (with 6" / 15.24 cm dual casters)	12.7" / 32.3 cm
High (with 6" / 15.24 cm dual casters)	31.7" / 80.5 cm
Patient surface positioning	
Backrest	0 ° to 60 °
Knee gatch	0 ° to 30 °
Trend. / Reverse Trend.	± 16 °
Boostless™	5" / 12.7 cm
Patient surface dimensions	
Width	35" to 47" / 89 cm to 119.4 cm
Length	80" or 84" / 203 cm or 213 cm
Extended Length	Adds 4" / 10.1 cm
Caster diameter	5" / 12.7 cm or 6" / 15.24 cm

Compliant with standards IEC 60601-1 3rd edition and IEC 60601-2-52. Certified for use by two people (referring to one patient and one bystander.)

Electronics

Volts	120 V AC
Ampere rating	4.25A (9.25A with 120V outlets)
Current leakage	< 300 micro amps
Frequency	60 Hz
Hospital grade plug	Yes
Battery type	Sealed rechargeable lead acid 1.2Ah, 2 * 12V
Scale system accuracy (in all bed positions)	1% with a maximum difference of 2.2 lb / 1 kg

Umano Medical reserves the right to change specifications without notice.

Umano Medical Inc.

230, boulevard Nilus-Leclerc, L'Islet QC G0R 2C0 CANADA

Toll free number : 1 844-409-4030

information@umanomedical.com

umanomedical.com