



umano medical

AN ALL-IN-ONE HOSPITAL BED FOR MULTIPLE CARE ENVIRONMENTS

med-surg | bariatric care | palliative care



 **ook snow** » ALL + **Pilot Drive**™

One-step Infinite™ width expansion system, increased weight capacity, and low height all come together to improve staff efficiency, patient comfort, and safety.

Featuring Pilot Drive™ – the first powered drive system to offer an unmatched combination of low height, underbed clearance, and maneuverability.

VERSATILITY • LOW HEIGHT • MOBILITY



MEDICAL-SURGICAL
Low height, ease of use, and ease of cleaning makes this bed ideal for med-surg environments.



BARIATRIC CARE
1060 lb weight capacity, width expansion, and an appealing design encourages patient autonomy and dignity for bariatric care.



PALLIATIVE CARE
Also known as the Cuddle Bed, it is the only bed on the market certified to be used by two people.



PILOT DRIVE
Added mobility combined with enhanced versatility makes the ook snow ALL with Pilot Drive and integrated O2 bottle holders a unique solution for your bariatric patients in critical care.



OPTIONAL FEATURES

- 11" low height with 5" dual casters
- Infinite™ width adjustment
- Integrated scale system
- Smart screen, nurse controls and patient control lockout on footboard
- Integrated bed exit detection system
 - 3 levels of sensitivity
 - Bed exit side view lights
- Smart Nightlights™
- Remote bed exit alarm connectivity
- Boostless™ Backrest System
- One hand release siderails
- Lifiable and flat deck covers
- Battery back-up
- Compatibility with mattresses within 48" width
- Programmable lowest height
- Nurse and patient controls integrated into siderails
- Integrated bed extender (35" x 4")
- CPR release handle
- Trendelenburg angle displayed on footboard
- Backrest angle indicator
- Patient restraint loops
- Integrated pump holder
- 120-volt / 5A insulated outlets (at foot end)
- Communication port (for bed software updates)
- Programmable Auto-Contour mode
- One-touch button groups
 - 30-degree head of bed
 - One-button chair position
 - Flat bed position
 - Trendelenburg / Reverse Trendelenburg
 - Vascular foot position

OPTIONAL FEATURES

- 12.7" low height 6" dual casters
- Pilot Drive™ with or without O2 bottle holders (compatible with size E bottles)
- Downward motion interrupt
- Bed status monitoring
- Optional patient surface length (80" or 84")
- 120-volt/5A insulated outlets with USB mobile device power outlets
- Higher siderails (for 8" mattress)
- Foot end back-up control on the frame
- Roller type wall protection bumpers
- Additional nurse call port
- One-touch button groups
 - Backlit patient controls
 - Nurse call
- Different styling options available

ACCESSORIES

- Patient helper system (300 lb capacity)
- Upright oxygen bottle holder
- IV poles (permanent & removable)
- 4" bed extender cushion
- Removable patient controls
- Transport cart

SPECIFICATIONS

ooksnow » ALL

Product model	200-0000
Length	
Overall length (80" deck length)	90" / 228.6 cm
Overall length (84" deck length)	94" / 239 cm
Extended length	Adds 4" / 10.1 cm
with Pilot Drive	Adds 1" / 2.5 cm
with Pilot Drive / oxygen bottle holders	Adds 2.25" / 5.7 cm
Overall width	
Siderails down - 35" surface	39.8" / 101 cm
Siderails up - 35" surface	41.5" / 105.4 cm
Siderails down - 47" surface	51.5" / 130.8 cm
Siderails up - 47" surface	52.8" / 134 cm
Safe Working Load	1060 lb / 481 kg
Height range (to deck top)	
Low (with 5" / 12.7 cm dual casters)	11" / 27.9 cm
High (with 5" / 12.7 cm dual casters)	30" / 76.2 cm
Low (with 6" / 15.2 cm dual casters)	12.7" / 32.3 cm
High (with 6" / 15.2 cm dual casters)	31.7" / 80.5 cm
Patient deck positioning	
Backrest	0° to 60°
Knee gatch	0° to 30°
Trend. / Reverse Trend.	± 16°
Boostless™	5" / 12,7 cm
Patient deck dimensions	
Width	35" to 47" / 89 cm to 119.4 cm
Length	80" or 84" / 203 cm or 213 cm
Extended Length	Adds 4" / 10.1 cm
Caster diameter	5" or 6" / 12.7 cm or 15.2 cm
with Pilot Drive	6" / 15.2 cm

Compliant with standards IEC 60601-1 and IEC 60601-2-52.
Certified for use by two people (referring to one patient and one bystander.)

ELECTRONICS

Volts	120 V AC
Ampere rating	4.25A (9.25A with 120V outlets)
Current leakage	Compliant with IEC 60601-1
Frequency	60 Hz
Hospital grade plug	Yes
Battery type	Sealed rechargeable lead acid 1.2Ah, 2 * 12V
Battery type (Pilot Drive)	Sealed rechargeable lead acid 16.5Ah, 2 * 12V
Scale system accuracy	1% with a maximum difference of 2.2 lb / 1 kg (in all bed positions)

Umano Medical reserves the right to change specifications without notice.

Umano Medical Inc.

230, boulevard Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll free number : 1 844-409-4030
information@umanomedical.com
umanomedical.com