



umano**medical**

AIRE-TURN™

POWERED THERAPEUTIC
SUPPORT SURFACE



Prevention and treatment of pressure injuries • Lateral rotation

KEY BENEFITS

The **Aire-Turn** is a powered support surface designed to address multiple needs:

- Assists in the prevention and treatment of pressure injuries.
- Provides 1 in 2 deep air cell alternating low pressure therapy.
- On demand low air loss therapy is administered directly underneath the patient, through the triple-layered top cover, and assists in regulating the microclimate to keep the skin cool and dry.
- Provides continuous lateral rotation therapy with adjustable cycle time and rotation angle of up to 45 degrees.*

*Depending on weight, position, and morphology of the patient.

Therapies

- Low Air Loss
- Alternating Pressure
- Lateral Rotation
- Pulsation / Massage



PRODUCT FEATURES

ANGLED AIR CELLS

Provides support and a more comfortable and subtle turn.

CONTROL UNIT

The setup button allows the caregiver to quickly set patient weight and height for automatic pressure settings, select desired therapies, and adjust cycle time, lateral rotation angle, duration, and sequence.

TOP COVER

Triple-layered top cover is low friction and low shear force producing, highly vapor permeable, and anti microbial. Available in nylon or 4-way stretch polyurethane material.



TECHNICAL SPECIFICATIONS

Model	K-6 Series	
Standard Size	36" W x 80" or 84" L x 10" H (91.4 cm W x 203 cm or 213.6 cm L x 25.4 cm H)	
Weight Capacity	500 lb (227 kg)	
Mattress Weight	31 lb (14 kg)	
Control Unit	Dimensions	12" x 6.75" x 10,5" (30 cm x 17 cm x 27 cm)
	Weight	13 lb (5.9 kg)
	Cycle Times	Alternating pressure cycle time: 1 - 99 minutes Lateral rotation cycle time: 10 minutes - 4 hours
	Power Cord	14' hospital grade detachable
	Input Voltage	90 ~ 240 VAC, 50/60 Hz, 1A/0.5A, Universal
Flammability	UL Listing	UL STD 60601-1 with respect to Electrical Shock, Fire and Mechanical Hazards
	Certified to	CAN/CSA STD C22.2 No. 601.1
Flammability	Mattress components meet flammability standard CAL 117 Optional Cal 129, 16 CFR 1632, CFR 1633	

Specifications can change without notice.



Distributed by Umano Medical Inc.

230, boulevard Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll Free: 1 (844) 409-4030
information@umanomedical.com
umanomedical.com

Manufactured by KAP Medical

1395 Pico Street Corona, CA USA 92881