



umano medical

THERAPEUTIC SERIES

airnest orange™



## Prevention and treatment of pressure injuries

### KEY BENEFITS

- The low air loss therapy can be activated or deactivated as needed.
- Offers real-time pressure monitoring sensor technology for accurate patient interface pressures and patient comfort.
- The system has lock out capabilities and internal memory that retains current settings.
- Holds air during transport and power outage.

### Therapies

- On-Demand Low Air Loss
- Alternating Pressure
- Static Mode



**PRODUCT FEATURES**

**DEEP AIR CELLS**

Pressure redistribution with alternating low-pressure therapy

**CONTROL UNIT**

Intuitive control panel allows the caregiver to quickly set patient weight, select desired therapies, and determine AP cycle times

**FOWLER BOOST MODE**

In upright position, the system increases the pressure in the air cells to prevent “bottoming out”



**TECHNICAL SPECIFICATIONS**

<b>Dimensions</b>	35" W x 80" or 84" L x 8" H (89 cm W x 203 cm or 213.4 cm L x 20.3 cm H)	
<b>Weight Capacity</b>	500 lb (227 Kg)	
<b>Mattress weight</b>	23 lb (10.4 Kg)	
<b>Control Unit</b>	<b>Dimensions</b>	13.8" x 11" x 6.5" (35 cm x 28 cm x 16.4 cm)
	<b>Weight</b>	10.2 lb (4.6 Kg)
	<b>Cycle Times</b>	AP Cycle Times: 5, 10, 15, & 20 minutes
	<b>Power Cord</b>	14' hospital-grade detachable
	<b>Input Voltage</b>	90 264 V, 50/60 Hz, 1A/0.5A, Universal
	<b>ETL Listed</b>	Conforms to UL STD 60601-1 with respect to Electrical Shock, Fire and Mechanical Hazards
	<b>Certified to</b>	CAN/CSA STD C22.2 No. 601.1
<b>Flammability</b>	Meets flammability standard CAL 117 Foam components meet CAL 129, 16 CFR 1632 & 16 CFR 1633	

Specifications can change without notice.



Manufactured for:  
**Umano Medical Inc.**  
 230, boul. Nilus-Leclerc  
 L'Islet QC G0R 2C0 CANADA  
 Toll Free: 1 (844) 409-4030  
 information@umanomedical.com  
 umano-medical.com