



umano medical



A LOW-HEIGHT HOSPITAL BED
FOR MEDICAL-SURGICAL AND HIGHER ACUITY CARE UNITS



Designed with a 360° approach on fall prevention and infection control, our low-height hospital bed will satisfy your most demanding needs in a way that's remarkably intuitive. Whether you're in a med-surg or high-acuity care unit, the **ooksnow** will redefine even your highest standards.

STANDARD BED FEATURES

- 10" low height
- Integrated scale
- Built-in bed exit detection system with 3 levels of sensitivity, Auto Arm, and Care Pause
- 5" dual casters
- Smart screen on footboard
- 3-position integrated bed extender
- Central braking system
- Bed exit side view lights
- Boostless™ Backrest System
- One hand release siderails
- Integrated patient and nurse controls with nurse call and backlight
- Nurse call configuration
- One-touch position buttons (Trendelenburg, vascular foot, chair, 30-degree head of bed and flat bed position)
- Backrest angle indicator
- Motion control master lockout
- 4 corner roller bumpers
- Integrated pump holder
- Roller type wall protection bumpers
- CPR release handle
- Quick-disconnect power cord
- Auxiliary outlets at foot end and side with USB ports
- Additional nurse call port
- Communication port (for easy bed software updates)
- Smart Nightlights™
- Battery back-up



Innovative hospital beds and support surfaces paired with unmatched service.

CHOOSE A BED PACKAGE ADAPTED TO YOUR CLINICAL NEEDS



PACKAGE UM-01 LOW HEIGHT AND SCALE

Designed with fall prevention in mind, this bed offers low height, a high-precision integrated scale and bed exit, without any compromise.

INCLUDED FEATURES

- ✓ Standard bed features
- ✓ 10" low height

UPGRADES AVAILABLE

- + Higher siderails (for 8" mattress)
- + Room lights on patient controls



PACKAGE UM-02 MONITORING

Designed for advanced local and remote monitoring of bed status to support your fall prevention strategy.

INCLUDED FEATURES

- ✓ Standard bed features
- ✓ 10" low height
- ✓ Bed status monitoring (brakes, siderails, backrest and bed height)

UPGRADES AVAILABLE

- + Umano Connect™ Technology (WiFi/BLE)
- + Higher siderails (for 8" mattress)
- + Room lights on patient controls



PACKAGE UM-03 MANEUVERABILITY

Designed to increase maneuverability without compromising low height and underbed clearance.

INCLUDED FEATURES

- ✓ Standard bed features
- ✓ Features of UM-02
- ✓ 12" low height
- ✓ 6" dual casters
- ✓ Pilot Drive™ powered drive system with two-handle control
- ✓ Integrated oxygen bottle holders

UPGRADES AVAILABLE

- + Umano Connect™ Technology (WiFi/BLE)
- + Higher siderails (for 8" mattress)
- + Room lights on patient controls

OPTIONAL FEATURES (CAN BE ADDED TO ANY BED PACKAGE)

- 2 or 5 function removable patient control
- Downward motion interrupt
- + Foot end back-up control on frame



CUSTOMIZE YOUR BED

The ook snow can also be built à la carte upon request according to the particular requirements of our customers. For more information, please contact your local account executive.

Discover our complete line of preventive and therapeutic **support surfaces** and all **accessories** available to enhance your bed functionalities on our website.

TECHNICAL SPECIFICATIONS

Product model	FL36	Patient deck positioning	0° to 60°
Weight capacity	600 lb / 272.5 kg	Backrest	0° to 30°
Caster diameter	5"	Knee gatch	± 16°
With Pilot Drive™	6"	Trendelenburg / Reverse Trend.	14° / 16°
		With Pilot Drive™	4"
		Boostless™ Backrest System	
Compliant with standard	IEC 60601-2-52	Electronics	IEC 60601-1 3rd edition
Patient deck dimensions	35" x 80" - 84" - 90"	Compliant with standard	120V AC
		Volts	4.75A (9.75A with 120V AC auxiliary outlets)
Overall length (80" deck)	88"	Current leakage	Conform with standards IEC 60601-1
Overall length extended (84" - 90" deck)	92" - 98"	Frequency	60 Hz
With roller type wall protection bumpers	Adds 7/8"	Hospital grade plug	Yes
With Pilot Drive™ / roller type wall protection bumpers	Adds 3.9"	Battery type	Sealed rechargeable lead acid battery 1.2 Ah, 2 * 12V
With Pilot Drive™ / integrated oxygen bottle holders / roller type wall protection bumpers	Adds 4.7"	Battery type (Pilot Drive™)	Sealed rechargeable lead acid battery 16.5 Ah, 2 * 12V
Overall width – siderails down	39.25"	Scale system accuracy	0.5 % with a maximum difference of 1.1 lb / 0.5 kg
Overall width – siderails up	41.3"	(in all bed positions)	
Height range (to deck top)			
Low (with 5" dual casters)	10"		
High (with 5" dual casters)	30"		
Low (with 6" dual casters)	12"		
High (with 6" dual casters)	31.8"		

Umano Medical reserves the right to change specifications without notice.



Umano Medical Inc.
230, boulevard Nilus-Leclerc
L'Islet QC G0R 2C0 CANADA
Toll free number : 1 844-409-4030
information@umanomedical.com
umanomedical.com

74-36964-ENG R2 USA - © Umano Medical Inc. 2023 - Printed in Canada

Watch the ook snow
in-service video

